

Manual Hydraulic Operating Table

MT300D

MT300D is widely used in chest, Abdominal surgery, ENT, gynaecology and obstetrics, urology and orthopedics etc.
Hydraulic lift by foot pedal, head operated movements.
The base and column cover all made by premium 304 stainless steel.



Model	Site	Length	Width	Tabletop lowest height	Elevation range
MT300D		≥1940mm	480mm	≤770mm	250mm

PARAMETER MT300D Manual Hydraulic Operating Table

Movements		Model:MT300D
Tabletop trendelenburg	Reverse Trendelenburg	≥25°
	Trendelenburg	≥25°
Lateral Tilt (Left And Right)		≥20°
Head plate	Up	≥60°
	Down	≥90°
Back plate	Up	≥80°
	Down	≥10°
Leg plate	Up	≥15°
	Down	≥90°
	Outward	≥90°
Kidney Bridge Elevation		≥120 mm